

EBV-VCA IgG Elisa

CAT NO	DESCRIPTION	PACK SIZE
EIAVGG1	EBV-VCA IgG Elisa	96 Tests
Intended Use:		
The Epstein-Barr Virus-Viral Capsid Antigen (EBV-VCA) IgG Enzyme linked Immunosorbent Assay (ELISA), is intended for the detection and quantitative determination of IgG antibody to Epstein-Barr virus in human sera. A single serum specimen may be used to indicate previous infection or immune status with the Epstein-Barr virus.		
Summary and Principle:		
Detection of the Epstein-Barr virus was first described in 1964 by Epstein, Achong, and Barr using electron microscopic studies of cultured lymphoblasts derived from patients with Burkitt's lymphoma. EBV is classified as a member of the herpes-virus family based upon its characteristic morphology. EBV infection may demonstrate a wide spectrum of clinical symptoms. The majority of primary EBV infections are transmitted via saliva, occur during childhood, and are subclinical. In the U.S., 50% of the population demonstrate EBV antibodies before the age of 5 years; 80% by adulthood. Transfusion-associated EBV infections have also been reported. In young adults, EBV infection may be clinically manifested as Infectious Mononucleosis (IM) with typical symptoms of sore throat, fever, and lymphadenopathy. College students and military personnel are often cited as a high morbidity incidence population for IM. Following primary EBV infection, it is postulated that the B lymphocyte may continue to harbor the EBV genome and establish a latent infection that may extend through life. Reactivation of EBV infection or enhanced EBV activation has been documented in immunodeficient or immunosuppressed patients, i.e., organ transplant patients, individuals with malignancies, pregnant women, and persons of advanced age. Epstein-Barr virus has also been associated in the pathogenesis of two human cancers, Burkitt's lymphoma and nasopharyngeal carcinoma. Documentation by means of DNA hybridization studies demonstrates the presence of the EBV genome on biopsy specimens taken from individuals with these carcinomas. Burkitt's lymphoma is primarily observed in Sub-Saharan Africa, especially in African children, and in New Guinea. Malarial infections are usually diagnosed in Burkitt's lymphoma patients and are suggested to be a co-factor. Nasopharyngeal carcinoma is observed in Asia, most notably in Southern China, and may have genetic or environmental influences as the co-factor. In the last two decades, serological methods have progressed from testing for the presence of non-specific heterophile antibodies to measuring levels of IgG or IgM formed against subunits of EBV antigen complexes. One of the best indicators of active EBV infection is antibody to viral capsid antigens, structural proteins necessary for replication of the virus. Viral capsid antigens are present in every cell infected with EBV. The IgM response to VCA is among the earliest detectable humoral immune responses, usually present at the onset of the disease and peaking within four to six weeks. VCA-IgM levels are also transient, declining rapidly and usually becoming undetectable within two to three months from onset of clinical symptoms.		
	Wash Concentrate (10x)	100ml
	TMB Substrate	12ml
	Enzyme conjugate	12ml
	Stop Solution	12ml
		PBS-Tween at pH 7.4. 10X concentrate. The concentrate must be diluted with 900 ml of distilled water before use. Once diluted it is stable at room temperature for two months.
		Amber bottle. Mixture of TMB and Hydrogen Peroxide solution. Ready to use. Once open, stable for 2 months at 2-8°C.
		Red colour solution. 1 vial containing 12ml of HRP labelled Anti Human IgG antibodies in Buffered saline. Once open, stable for 2 months at 2-8°C.
		Diluted Sulfuric acid solution (1M) Ready to use. Once open, stable for 2 months at 2-8°C.

Plastic Sealable bag, IFU and Cardboard plate covers.

Specimen Collection:

1. Collect blood specimens and separate the serum.
2. Specimens may be refrigerated at 2 - 8° C for up to seven days or frozen for up to six months. Avoid repetitive freezing and thawing of serum sample.
3. If paired sera are to be collected, acute samples should be collected as soon as possible after onset of symptoms and not later than seven days after onset. The second sample should be collected 14 to 21 days after the acute specimen was collected. Both samples must be run in duplicate tests on the same plate to test for a significant rise. If the first specimen is obtained too late during the course of the infection, a significant rise may not be detectable.

Storage and Stability:

1. Store the kit at 2 - 8°C.
2. Always keep microwells tightly sealed in pouch with desiccants. We recommend you use up all wells within 4 weeks after initial opening of the pouch.
3. The reagents are stable until expiration of the kit.
4. Do not expose test reagents to heat, sun, or strong light during storage or usage.

Precautions and Safety:

1. Potential biohazardous materials: The calibrator and controls contain human source components, which have been tested and found nonreactive for Hepatitis B surface antigen as well as HIV antibody with FDA licensed reagents. However, as there is no test method that can offer complete assurance that HIV, Hepatitis B virus, or other infectious agents are absent, these reagents should be handled at the Biosafety Level 2, as recommended in the Centers for Disease Control / National Institutes of Health manual, "Biosafety in Microbiological and Biomedical Laboratories." 1984
2. Do not pipette by mouth. Do not smoke, eat, or drink in the areas in which specimens or kit reagents are handled.
3. The components in this kit are intended for use as an integral unit. The components of different lots should not be mixed.
4. This product contains components preserved with sodium azide. Sodium azide may react with lead and copper plumbing to form explosive metal azide. On disposal, flush with a large volume of water.

Procedure:

Reagent preparation:

Prepare 1x washing buffer. Prepare washing buffer by adding distilled or deionized water to 10x wash concentrate to make a final volume of 1 litre.
Bring all specimens and kit reagents to room temperature 20-25°C and gently mix.

STEP 1

Preparation: Place the desired number of coated strips into the holder. Pre-wash Coated wells – repeat washing three times with washing buffer.

STEP 2

Addition of the Diluent: Prepare 1:21 dilutions by adding 10ul of the samples, negative control, positive control, and calibrator to 200 ul of sample diluent. Mix well.

STEP 3

Addition of the Sera, calibrators and controls: Dispense 100ul of diluted sera, calibrator and controls into the appropriate wells. For the reagent blank, dispense 100 ul sample diluent in 1A well position. Tap the holder to remove air bubbles from the liquid and mix well.

STEP 4

Incubation: Incubate for 30 minutes at room temperature.

STEP 5

COMPONENT	SIZE	DESCRIPTION
Microwell Plate	1x96 wells (12x8 well plate)	Each microwell is coated with Purified EBV VCA Antigens. The microwells can be broken and used separately. Place unused wells or strips in the provided plastic sealable bag together with the desiccant and store at 2-8°C. Once open the wells are stable for 2 months at 2-8°C.
EBV VCA IgG Calibrator	150ul	Factor value (f) stated on label. Red cap. Ready to use. Once open, stable for 1 month at 2-8°C.
Negative Control	150ul	Range stated on label. Natural cap. Control material negative for EBV VCA IgG. Ready to use. Once open, stable for 1 month at 2-8°C.
Positive Control	150ul	Range stated on label. Green cap. Control Positive for EBV VCA IgG. Ready to use. Once open, stable for 1 month at 2-8°C.
Sample Diluent	22ml	Blue colour solution. 1 vial containing 22 ml of Sample Diluent made from Buffer. Store at 2-8°C. Once opened stable for 2 months at 2-8°C.

Washing: Remove liquid from all wells. Repeat washing three times with washing buffer.

STEP 6

Addition of Enzyme Conjugate: Dispense 100ul of enzyme conjugate to each well.

STEP 7

Incubation: Incubate for 30 minutes at room temperature.

STEP 8

Washing: Remove enzyme conjugate from all wells. Repeat washing three times with washing buffer.

STEP 9

Addition of TMB Chromogenic Substrate: Dispense 100ul of TMB Substrate into each well and incubate for 15 minutes at room temperature.

STEP 10

Addition of Stop solution: Add 100ul of Stop Solution to stop reaction. Make sure there are no air bubbles in each well before reading.

STEP 11

Measurement: Read O.D. at 450nm with a microwell reader.

Calculation of results:

- To obtain Cut off OD value: Multiply the OD of Calibrator by Factor (f) printed on label of Calibrator.
- Calculate the EBV-VCA IgG Index of each determination by dividing the OD values of each sample by obtained OD value of Cut off.

For example:

- If Factor (f) value on label = 0.4
- This factor (f) is a variable. It is specific for a lot manufactured and printed on label of Calibrator.
- Obtained Calibrator O.D. = 1.100.
- Cut-off O.D. = $1.100 \times 0.4 = 0.44$ (By definition EBV-VCA IgG Index = 1)
- Patient sample O.D. = 0.580
- EBV-VCA IgG Index = $0.580 / 0.44 = 1.32$ (Positive result)
- Patient sample O.D. = 0.320
- EBV-VCA IgG Index = $0.320 / 0.44 = 0.73$ (Negative result)

Quality Control:

The test run may be considered valid provided the following criteria are met:

- The O.D. value of the reagent blank against air from a microwell reader should be less than 0.150.
- If the O.D. value of the Calibrator is lower than 0.30, the test is not valid and must be repeated.
- The EBV-VCA IgG Index for Negative and Positive Control should be in the range stated on the labels.
-

Interpretation of results:

Negative: EBV-VCA Index of 0.90 or less are seronegative for IgG antibody to EBV-VCA virus.

Equivocal: EBV-VCA Index of 0.91 - 0.99 are equivocal. Sample should be retested.

Positive: EBV-VCA Index of 1.00 or greater are seropositive. It indicates prior exposure to the EBV-VCA virus.

Significant change in antibody level of the paired sample:

The ratio between the EBV-VCA IgG Index of the second sample and that of the first sample should be greater than 1.3 to be suggestive of a significant rise in antibody level.

Limitation of the Test:

- A single serum sample cannot be used to determine recent infection.
- A serum specimen taken in an early stage during acute phase of infection may contain low levels of IgG antibody and render an EBV-VCA IgG Index result negative.
- As with other serological assays, the results of these assays should be used in conjunction with information available from clinical evaluation and other diagnostic procedures.

Performance Characteristics:

Sensitivity, specificity and accuracy:

A total of 225 random samples from different sources were assayed with EBV-VCA IgG test and with another commercially available Elisa test kit.


REFERENCE ELISA					
		N	E	P	TOTAL
	N	11 (D)	0	0 (B)	11
Prestige Elisa	E	0	3	0	3
	P	0 (C)	0	211 (A)	211
	TOTAL	11	3	211	225

Sensitivity = $A / (A+B) = 211 / 211 = 100\%$

Specificity = $D / (C+D) = 11 / 11 = 100\%$

Accuracy = $(A+D) / (A+B+C+D) = 222 / 222 = 100\%$

The correlation of quantitative values between two comparison methods was summarized:



Precision:

The mean, SD and % CV were calculated inter- and intra- assay:

Intra-assay	n	Mean	SD	%CV
Serum 1	8	0.1705	0.0074	4.39
Serum 2	8	1.274	0.0997	7.83
Serum 3	8	2.330	0.1640	7.04
Serum 4	8	2.481	0.1099	4.43
Inter-assay	n	Mean	SD	%CV
Serum 1	8	0.170	0.0071	4.20
Serum 2	8	1.296	0.0486	3.75
Serum 3	8	2.278	0.1220	5.36
Serum 4	8	2.444	0.2543	10.41

Limitations of the Assay:

- The antibody titer of a single serum specimen cannot be used to determine recent infection. Test results for anti-VCA should be interpreted in conjunction with the clinical evaluation and results of antibody tests for other EBV antigens.
- Most (80%) of IM individuals have peak anti-VCA IgG titers before they consult a physician. Therefore, testing paired acute and convalescent sera for significant changes in antibody levels is not useful in most patients with IM.

References:

- Epstein, M.A., B.B. Achong, and Y.M. Barr. 1964. Virus particles in Cultured Lymphoblasts from Burkitt's Lymphoma. In: Lancet 1:702- 703.
- Epstein, M.A., Y.M. Barr, and B.G. Achong. 1965. Studies with Burkitt's Lymphoma. In: Wistar Inst. Sympos. Monogr. 4:69-82.
- Schooley, R.T. and R. Dolin. 1985. Epstein-Barr Virus Infectious Mononucleosis. In: Principles and Practice of Infectious Diseases, 2nd Edition. Mandell, G.L., R.G. Douglas, and J.E. Bennett. (eds). John Wiley and Sons, New York. pp 97-982.
- Sumaya, C.V. 1985. Serological Testing for Epstein-Barr Virus- Developments in Interpretations. In: J. Inf. Dis. 151 (6):984-987.
- Tobi, M., and S.E. Straus. 1985. Chronic Epstein-Barr Virus Disease: A Workshop Held by the National Institute of Allergy and Infectious Diseases. In: Ann. Intern. Med. 103 (6 (pt. 1)):951-953.
- Birx, D.L., R.R. Redfield, and G. Tosarto. 1986. Defective Regulation of Epstein-Barr Virus Infection in Patients with Acquired Immunodeficiency Syndrome (AIDS) or AIDS-Related Disorders. In: New Eng. J. Med. 314 (14):874-879.
- Pearson, G.R. and Luka, J. Characterisation of the Virusdetermined Antigens. In The Epstein-Barr Virus: Recent Advances. Epstein and Achong, eds. Wiley: New York. 1986.
- Schmitz, H. and Scherer, M. IgM Antibodies to Epstein-Barr Virus in Infectious Mononucleosis. Arch. Ges. Virusforsch. 37:332-39, 1972.

REF	Catalog number	Lot	Temperature limitation
	Consult instructions for use		Batch code
	In vitro diagnostic medical device		Use by
	Manufacturer		

