


Dengue Dual Ag/Ab Rapid device

Innovative life solutions from Prestige



Primary Infection Pathway of the Dengue Virus

Dengue Dual Ag/IgG-IgM Rapid Device

For testing Serum / Plasma / Whole Blood

- ◆ 20 minutes testing time
- ◆ Detects all 4 types of Dengue serotype.
- ◆ For Detection of current infection and recent infection
- ◆ Detects as rapidly as 2 days after infection.
- ◆ Excellent Clinical performance
- ◆ Antigen Sensitivity: 96.4%, Antigen Specificity: 96.0%
- ◆ IgM Sensitivity 96.0%, IgM Specificity: 97.3%
- ◆ IgG Sensitivity 95.6%, IgG Specificity: 96.1%
- ◆ Detects in all 3 phases of Dengue infection—Febrile Phase
Critical Phase and Recovery Phase.

Dengue infection

400 million infected every year


- ◆ 4 Serotypes of Dengue Virus
- ◆ DEN1, DEN2, DEN3, DEN4
- ◆ Diagnosis Goal: Early Detection
- ◆ There is no known treatment for Dengue.

Prestige Diagnostics Dengue Dual Device

*For early detection of
Dengue Fever.*

Prestige Diagnostics— House of Rapid Tests

Dengue, countries or areas at risk, 2011



Ordering Information: RADDEN1: 25 Tests

Prestige Diagnostics UK Ltd., Ballymena Business Centre, Galgorm, Co.Antrim, BT42 1FL, UK
+44-(0)28 25642100 www.prestigediagnostics.co.uk Email: info@prestigediagnostics.co.uk